

## Laser Ablation - Inductively Coupled Plasma - Mass Spectrometry (LA-ICP-MS)

### Technology:

**Mass Spectrometry of Material Ablated by Laser, Surface and Micro Area Analysis**

### Equipment:

- ICP-MS: 7500ce, Agilent
- UP 193 FX, New Wave

### Category:

**C. Particle Characterisation ex-situ**

### Institute:

**KIT**

### Location:

Karlsruhe Institute of Technology  
 Institute for Applied Materials (IAM-AWP)  
 Building 688  
 Hermann-v.-Helmholtz-Platz 1  
 D-76344 Eggenstein-Leopoldshafen  
 Germany

### Contact Details of Technology Expert:

**Name:** Dr. Thomas Bergfeldt

**Phone:** +49 721 608-22914

**Fax:** +49 721 608-922914

**E-mail:** Thomas.Bergfeldt@kit.edu

### Short technology description/Overview and Main Features (Equipment Capabilities)::

LA-ICP-MS is a sensitive analytical method (0,000001 – 100%) for rapid multi-element (Li - U except C, N, O, F, Cl) determination in the trace and ultra trace range of different solid sample materials, technical products and even biological samples (biological samples can't be analyzed in this laboratory). A small area ( $\varnothing$  10 – 200  $\mu$ m) of the sample is vaporized in a laser plasma by focused laser radiation and transported with helium in the inductively coupled plasma ion source of an ICP-MS. There the material is atomized, ionized, accelerated into the mass spectrometer and separated according the mass/charge ratio and energy/charge ratio and detected by electron multipliers.

The main advantage of LA-ICP-MS is that samples are investigated not under vacuum but in atmospheric pressure and electrically conducting as well as non-conducting material can be investigated.

### Typical Samples & results:

Quantitative analysis of small metal fragments by LA-ICP-MS. Instrument calibration was achieved using 8 flat geometry steel setting-up samples (BAS UK). The samples included low-alloy and highly-alloyed steels. The elements varied from ppm to high % levels.

	Al		P		V		Cr		Nb		Mo		Sn	
	Cert	LA	Cert	LA	Cert	LA	Cert	LA	Cert	LA	Cert	LA	Cert	LA
A/8	.04	.04	.005	.01	.005	.001	.02	.03	.01	---	.005	.002	.005	.001
D/11	.19	.19	.01	.02	.12	.11	3.0	3.1	.05	.05	1.3	1.3	.01	.01
C/17	.05	.05	.07	.08	.42	.42	.18	.20	.02	.03	.11	.11	.06	.07
H/5	.2	.16	.04	.06	.33	.33	1.3	1.1	.11	.09	.42	.45	.03	.03
G/5	.005	.006	.03	.04	.06	.07	16.5	12.9	.005	.003	2.2	2.2	.01	.02
E/4	.03	.04	.01	.02	.03	.04	4.3	3.6	1.1	1.2	.97	.98	---	.002
B/6	.01	.02	.005	.01	.01	.01	1.2	1.2	.005	.006	.20	.21	.01	.01
F/4	.005	.004	.02	.03	.05	.05	25.5	23.4	.005	.001	3.5	3.4	.01	.02

Table 1: Comparison of LA-ICP-MS values with certified values (bulk samples). All data % (w/w)

Acknowledgements to Dr. Walter Grillo and Dr. Massimo Godano, FIAT, Turin, Italy

Any further Information: